

VisTrill[®] INJECTOR AND CARTRIDGE

HEALTH CANADA APPROVED DEVICE ID NUMBER : 1036232

TMC 503 flat-back loaded cartridge technology allows a surgical team to smoothly manoeuvre an IOL into the cartridge with ease. When the IOL is safely manipulated into place, the cartridge's lock design ensures sleek and secure loading for safe lens delivery.

- Perfect for mono manual and bi-manual injection procedures
- Flat-back loading cartridge for easy lens loading and safe lens manipulation for the end user
- Safe lens delivery with a cartridge locking system



DELIVERY SYSTEM
AVAILABLE IN 2.2 mm AND 2.4 mm



Wound-assisted technique
External Diameter: 1.6 - 2.4



Into-wound-assisted technique
External Diameter: 1.8 - 2.8



Into-bag-assisted technique
External Diameter: 2.0 - 3.0

The incision size depends on the ophthalmic surgeon technique.